


Anti-Traf3 Rabbit Polyclonal Antibody For Wb - Q13114

Item Number: CM0029221



Introduction

TRAF3 rabbit polyclonal antibody (pAb) for WB, IHC-P, IF/ICC, ELISA. Reacts with Human, Mouse, Rat. Validated for cytoplasmic E3 ubiquitin ligase research in immune signaling.

[Learn More](#)

Parameter	Value
Product Name	TRAF3 Rabbit pAb
Remarks/Alias	CAP1; LAP1; CAP-1; CRAF1; IIAE5; CD40bp; RNF118; TRAF3
Species	Human
Gene ID (Human)	7187
Gene ID	7187
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human TRAF3 (NP_003291.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13114
Protein Weight	64kDa
Shipping	Ice bag