

Anti-Traf3Ip2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O43734

Item Number: CM0030992



Introduction

High-quality rabbit polyclonal antibody to TRAF3IP2 (ACT1/CIKS), an E3 ubiquitin ligase essential for IL-17 signaling and NF-κB activation. Suitable for WB, IF/ICC, ELISA in human, mouse, rat.

[Learn More](#)

Parameter	Value
Product Name	TRAF3IP2 Rabbit pAb
Remarks/Alias	ACT1; CIKS; C6orf2; C6orf4; C6orf5; C6orf6; CANDF8; PSORS13; TRAF3IP2
Species	Human
Gene ID (Human)	10758
Gene ID	10758
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human TRAF3IP2 (NP_671733.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43734
Protein Weight	65kDa
Shipping	Ice bag