


Anti-Traf4 Rabbit Monoclonal Antibody For Wb, Elisa - Q9Buz4

Item Number: CM0006690



Introduction

Rabbit monoclonal antibody against human TRAF4, validated for WB and ELISA, cross-reacts with mouse and rat. Suitable for studying cell proliferation, differentiation, DNA repair, and NF-κB signaling.

[Learn More](#)

Parameter	Value
Product Name	TRAF4 Rabbit mAb
Remarks/Alias	CART1; MLN62; RNF83; TRAF4
Species	Human
Gene ID (Human)	9618
Gene ID	9618
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 209-470 of human TRAF4 (NP_004286.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BUZ4
Protein Weight	54 kDa
Shipping	Ice bag