


Anti-Traf6 Rabbit Monoclonal Antibody For Wb, Ip, Elisa - Q9Y4K3

Item Number: CM0007997



Introduction

High-quality anti-TRAF6 rabbit monoclonal antibody validated for WB, IP, ELISA in human, mouse, and rat. Targets TNF receptor-associated factor 6 (Q9Y4K3), a key E3 ubiquitin ligase in NF-κB signaling.

[Learn More](#)

Parameter	Value
Product Name	TRAF6 Rabbit mAb
Remarks/Alias	RNF85; MGC:3310
Species	Human
GeneID (Human)	7189
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-274 of human TRAF6 (NP_004611.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y4K3
Protein Weight	60 kDa
Shipping	Ice bag