


Anti-Traf6 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Y4K3

Item Number: CM0021268



Introduction

High-quality TRAF6 Rabbit Polyclonal Antibody, validated for WB, IF/ICC, ELISA in human, mouse, rat. Ideal for studying ubiquitin ligase function and NF-κB signaling pathways. Specific and reliable for research use.

[Learn More](#)

Parameter	Value
Product Name	TRAF6 Rabbit pAb
Remarks/Alias	RNF85; MGC:3310; TRAF6
Species	Human
Gene ID	7189
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 50-150 of human TRAF6 (NP_004611.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y4K3
Protein Weight	60kDa
Shipping	Ice bag