


Anti-Traf6Bp/Tax1Bp1 Monoclonal Antibody For Wb, Elisa - Q86Vp1

Item Number: CM0003477



Introduction

Recombinant rabbit monoclonal antibody targeting TRAF6BP/TAX1BP1 (T6BP/TXBP151). Validated for WB and ELISA with cross-reactivity to human and mouse. Ideal for innate immunity and autophagy studies.

[Learn More](#)

Parameter	Value
Product Name	TRAF6BP/TAX1BP1 Rabbit mAb
Remarks/Alias	T6BP; TXBP151; CALCOCO3; TRAF6BP/TAX1BP1
Species	Human
GeneID (Human)	8887
GeneID	8887
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 690-789 of human TRAF6BP/TAX1BP1 (Q86VP1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q86VP1
Protein Weight	91kDa
Shipping	Ice bag