


Anti-Traf6Bp/Tax1Bp1 Polyclonal Antibody For Ihc-P And Elisa - Q86Vp1

Item Number: CM0022950



Introduction

Rabbit polyclonal antibody targeting TRAF6BP/TAX1BP1 (T6BP, TXBP151, CALCOCO3). Suitable for IHC-P and ELISA, with cross-reactivity for human and rat. Ideal for research on NF-kB signaling, autophagy, and innate immunity.

[Learn More](#)

Parameter	Value
Product Name	TRAF6BP/TAX1BP1 Rabbit pAb
Remarks/Alias	T6BP; TXBP151; CALCOCO3; TRAF6BP/TAX1BP1
Species	Human
GeneID (Human)	8887
GeneID	8887
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 709-789 of human TRAF6BP/TAX1BP1 (NP_006015.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	Q86VP1
Protein Weight	91kDa
Shipping	Ice bag