

## Anti-Traf7 Polyclonal Antibody For Wb, Elisa - Q6Q0C0

Item Number: CM0013999



### Introduction

High-quality rabbit polyclonal antibody targeting human TRAF7 (Q6Q0C0) for WB and ELISA applications. Validated cross-reactivity: Human. Ideal for research in innate immunity, inflammation, and apoptosis pathways.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | TRAF7 Rabbit pAb                                                                                                   |
| Remarks/Alias    | RFWD1; RNF119; CAFDADD; TRAF7                                                                                      |
| Species          | Human                                                                                                              |
| Gene ID (Human)  | 84231                                                                                                              |
| Gene ID          | 84231                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 260-400 of human TRAF7 (NP_115647.2) |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, ELISA                                                                                                          |
| Cross Reactivity | Human                                                                                                              |
| SWISS            | Q6Q0C0                                                                                                             |
| Protein Weight   | 75kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                            |