

Anti-Traf7 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q6Q0C0

Item Number: CM0031023



Introduction

TRAF7 Rabbit pAb is a rabbit polyclonal antibody for WB, IHC-P, IF/ICC, ELISA applications. It reacts with human, mouse, rat TRAF7. Ideal for ubiquitin ligase, innate immunity, and apoptosis studies.

[Learn More](#)

Parameter	Value
Product Name	TRAF7 Rabbit pAb
Remarks/Alias	RFWD1, RNF119, CAFDADD, TRAF7
Species	Human
GeneID (Human)	84231
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6Q0C0
Protein Weight	75kDa
Shipping	Ice bag