

Anti-Traf7 Rabbit Polyclonal Antibody For Western Blot And Elisa - Q6Q0C0

Item Number: CM0024362



Introduction

TRAF7 rabbit polyclonal antibody validated for Western blot and ELISA. Detects human TRAF7 with cross-reactivity to rat. Ideal for investigating innate immunity, inflammation, and apoptosis pathways. Supplied by CamelBio.

[Learn More](#)

Parameter	Value
Product Name	TRAF7 Rabbit pAb
Remarks/Alias	RFWD1; RNF119; CAFDADD; TRAF7
Species	Human
Gene ID (Human)	84231
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 260-400 of human TRAF7 (NP_115647.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q6Q0C0
Protein Weight	75kDa
Shipping	Ice bag