


Anti-Trafd1 Rabbit Polyclonal Antibody For Wb, Elisa - O14545

Item Number: CM0013168



Introduction

TRAFD1 Rabbit Polyclonal Antibody detects human TRAFD1 (FLN29) via WB and ELISA. Target is a negative regulator of innate immune TLR4/RLH pathways, useful for inflammation and signaling research.

[Learn More](#)

Specification	Detail
Product Name	TRAFD1 Rabbit pAb
Remarks/Alias	FLN29; TRAFD1
Species	Human
Gene ID (Human)	10906
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human TRAFD1 (NP_006691.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O14545
Protein Weight	65kDa
Shipping	Ice bag