

Anti-Tram1 Rabbit Polyclonal Antibody For Wb - Q15629

Item Number: CM0013531

Introduction

Rabbit polyclonal antibody against human TRAM1, validated for WB and ELISA. Cross-reacts with mouse. Useful for studying ER protein translocation, unfolded protein response, and viral immune evasion.

[Learn More](#)

Parameter	Value
Product Name	TRAM1 Rabbit pAb
Remarks/Alias	TRAM; PNAS8; TRAMP; TRAM1
Species	Human
GeneID (human)	23471
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 305-374 of human TRAM1 (NP_055109.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q15629
Protein Weight	43kDa
Shipping	Ice bag