


Anti-Transcobalamin-1 (Tcn1) Polyclonal Antibody For Wb, Elisa - P20061

Item Number: CM0028443



Introduction

High-quality rabbit polyclonal antibody targeting human Transcobalamin-1 (TCN1/Haptocorrin). Validated for WB and ELISA, with cross-reactivity to mouse and rat. Ideal for studying vitamin B12 transport and gastrointestinal protection.

[Learn More](#)

Parameter	Value
Product Name	TCN1 Rabbit pAb
Remarks/Alias	HC, TC1, TCI, TC-1, TCN1
Species	Human
Gene ID (Human)	6947
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-280 of human TCN1 (NP_001053.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS Accession	P20061
Protein Weight	48kDa
Shipping	Ice bag