

Anti-Transferrin Polyclonal Antibody For Wb, If, Elisa - P02787

Item Number: CM0014425



Introduction

This rabbit polyclonal antibody targets human Transferrin (Sero-transferrin), validated for WB, IF/ICC, IF-P, and ELISA. It cross-reacts with mouse and rat, making it ideal for iron transport, cellular iron uptake, and homeostasis studies.

[Learn More](#)

Parameter	Value
Product Name	Transferrin Rabbit pAb
Remarks/Alias	TFQTL1; PRO1557; PRO2086; HEL-S-71p; Transferrin
Species	Human
Gene ID (Human)	7018
Gene ID	7018
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 355-698 of human Transferrin (NP_001054.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P02787
Protein Weight	77kDa
Shipping	Ice bag