

Anti-Transferrin Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P02787

Item Number: CM0002746



Introduction

This high-quality rabbit monoclonal antibody targets human Transferrin, validated for Western blot, IHC-P, and ELISA. It is an essential tool for iron metabolism, cancer, and infectious disease research.

[Learn More](#)

Parameter	Value
Product Name	Transferrin Rabbit mAb
Remarks/Alias	TFQTL1; PRO1557; PRO2086; HEL-S-71p; Transferrin
Species	Human
GeneID (Human)	7018
GeneID	7018
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 50-150 of human Transferrin (P02787)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P02787
Protein Weight	77kDa
Shipping	Ice bag