


Anti-Transglutaminase 2 (Tgm2) Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P21980

Item Number: CM0005489



Introduction

Rabbit monoclonal antibody against human Transglutaminase 2 (TGM2) validated for WB, IHC-P, IF/ICC, ELISA. Ideal for studying protein cross-linking, apoptosis, and chromatin modification. Research use only.

[Learn More](#)

Parameter	Value
Product Name	Transglutaminase 2 (TGM2) Rabbit mAb
Remarks/Alias	TGC; tTG; G(h); hTG2; TG(C); Transglutaminase 2 (TGM2)
Species	Human
GeneID (Human)	7052
GeneID	7052
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 438-687 of human Transglutaminase 2 (TGM2) (NP_004604.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P21980
Protein Weight	77kDa
Shipping	Ice bag