


Anti-Transglutaminase 2 (Tgm2) Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P21980

Item Number: CM0017608



Introduction

Transglutaminase 2 (TGM2) rabbit polyclonal antibody for Western blot, IHC-P, IF/ICC, and ELISA. Cross-reacts with human, mouse, rat. Used in apoptosis, chromatin modification, and extracellular matrix studies.

[Learn More](#)

Parameter	Value
Product Name	Transglutaminase 2 (TGM2) Rabbit pAb
Remarks/Alias	TGC; tTG; G(h); hTG2; TG(C); Transglutaminase 2 (TGM2)
Species	Human
GeneID (Human)	7052
GeneID	7052
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 299-548 of human Transglutaminase 2 (TGM2) (NP_945189.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P21980
Protein Weight	77kDa
Shipping	Ice bag