

Anti-Transketolase (Tkt) Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P29401

Item Number: CM0030874



Introduction

CamelBio offers a rabbit polyclonal antibody against human Transketolase (TKT), validated for WB, IHC-P, and ELISA with cross-reactivity to mouse and rat. Ideal for studying the pentose phosphate pathway and metabolic disorders.

[Learn More](#)

Parameter	Specification
Product Name	TKT Rabbit pAb
Remarks/Alias	TK; TKT1; SDDHD; HEL107; HEL-S-48; TKT
Species	Human
Gene ID (Human)	7086
Gene ID	7086
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human TKT (NP_001128527.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P29401
Protein Weight	68kDa
Shipping	Ice bag