

Anti-Translin (Tsn) Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q15631

Item Number: CM0023167



Introduction

Rabbit polyclonal antibody against human Translin (TSN, C3PO); validated in WB, IHC-P, IF/ICC, ELISA; cross-reacts with mouse and rat; detects ~26 kDa protein. Ideal for chromosomal translocation and RISC research.

[Learn More](#)

Parameter	Value
Product Name	TSN Rabbit pAb
Remarks/Alias	C3PO; RCHF1; TBRBP; TRSLN; BCLF-1; REHF-1; TSN
Species	Human
GeneID (Human)	7247
GeneID	7247
Immunogen	Synthetic peptide (a synthetic peptide corresponding to a sequence within amino acids 50-150 of human TSN, NP_004613.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15631
Protein Weight	26kDa
Shipping	Ice bag