

Anti-Transportin 3 (Tnpo3) Rabbit Monoclonal Antibody For Wb, If/Icc, If-P, Ihc-P, Elisa - Q9Y5L0

Item Number: CM0002275



Introduction

Rabbit monoclonal antibody against human Transportin 3 (TNPO3), validated for WB, IF/ICC, IF-P, IHC-P, and ELISA; reacts with human, mouse, and rat. Ideal for studying nuclear transport and HIV-1 infection mechanisms.

[Learn More](#)

Parameter	Value
Product Name	Transportin 3 (TNPO3) Rabbit mAb
Remarks/Alias	IPO12, TRNSR, LGMD1F, LGMD22, MTR10A, TRN-SR, TRN-SR2, Transportin 3 (TNPO3)
Species	Human
GeneID (Human)	23534
GeneID	23534
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 824-923 of human Transportin 3 (Transportin 3 (TNPO3)) (Q9Y5L0).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y5L0
Protein Weight	104kDa
Shipping	Ice bag