


Anti-Transthyretin (Ttr) Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P02766

Item Number: CM0015194



Introduction

Rabbit polyclonal antibody targeting human Transthyretin (TTR), validated for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for studies on thyroid hormone transport and amyloidosis research.

[Learn More](#)

Parameter	Value
Product Name	TTR Rabbit pAb
Remarks/Alias	CTS; TTN; ATTR; CTS1; PALB; TBPA; HEL111; HsT2651; TTR
Species	Human
GeneID (Human)	7276
GeneID	7276
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-147 of human TTR (NP_000362.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P02766
Protein Weight	16kDa
Shipping	Ice bag