

Anti-Trap1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q12931

Item Number: CM0025890



Introduction

TRAP1 Rabbit pAb is a polyclonal antibody targeting human TRAP1 (HSP75, HSP90L). Validated for WB, IHC-P, IF/ICC, ELISA with human, mouse, rat reactivity. Ideal for mitochondrial chaperone and cancer research.

[Learn More](#)

Parameter	Value
Product Name	TRAP1 Rabbit pAb
Remarks/Alias	HSP75; HSP 75; HSP90L; TRAP-1; TRAP1
Species	Human
Gene ID	10131
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q12931
Protein Weight	80kDa
Shipping	Ice bag