

Anti-Trap1 Rabbit Monoclonal Antibody For Wb And Elisa - Q12931

Item Number: CM0011147



Introduction

This rabbit monoclonal antibody specifically recognizes human TRAP1 (HSP75), validated for WB and ELISA. Reacts with human and mouse samples, making it suitable for mitochondrial chaperone studies and cancer research.

[Learn More](#)

Parameter	Value
Product Name	TRAP1 Rabbit mAb
Remarks/Alias	HSP75, HSP 75, HSP90L, TRAP-1, TRAP1
Species	Human
GeneID	10131
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human TRAP1 (Q12931)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q12931
Protein Weight	80kDa
Shipping	Ice bag