

Anti-Trap1 Rabbit Monoclonal Antibody For Wb, Ip, If/Icc, If-P, Ihc-P, Elisa - Q12931

Item Number: CM0012712



Introduction

Validated anti-TRAP1 rabbit monoclonal antibody, detects human and mouse HSP75. Suitable for WB, IP, IF, IHC-P, ELISA. Ideal for mitochondrial function and cancer metabolism research.

[Learn More](#)

Parameter	Value
Product Name	TRAP1 Rabbit mAb
Remarks/Alias	HSP75; HSP 75; HSP90L; TRAP-1
Species	Human
Gene ID (Human)	10131
Gene ID	10131
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 445-704 of human TRAP1 (NP_057376.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q12931
Protein Weight	80kDa
Shipping	Ice bag