

Anti-Trap1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc - Q12931

Item Number: CM0029452



Introduction

TRAP1 Rabbit pAb is a polyclonal antibody validated for WB, IHC-P, IF/ICC, and ELISA in human, mouse, rat. It targets the mitochondrial chaperone TRAP1, involved in mitochondrial function and cancer metabolism. Ideal for oncology and cell stress research.

[Learn More](#)

Parameter	Value
Product Name	TRAP1 Rabbit pAb
Remarks/Alias	HSP75; HSP 75; HSP90L; TRAP-1; TRAP1
Species	Human
Gene ID (Human)	10131
Gene ID	10131
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 445-704 of human TRAP1 (NP_057376.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q12931
Protein Weight	80kDa
Shipping	Ice bag