

Anti-Trappc2 Monoclonal Antibody For Wb, Ihc-P, Elisa - P0Di81

Item Number: CM0003156

Introduction

Rabbit monoclonal antibody targeting human TRAPPC2/Sedlin. Validated for western blot, IHC-P, and ELISA. Cross-reacts with mouse and rat. Suitable for vesicular transport and transcription studies.

[Learn More](#)

Parameter	Value
Product Name	TRAPPC2 Rabbit mAb
Remarks/Alias	SEDL; SEDT; MIP2A; TRS20; ZNF547L; hYP38334; TRAPPC2P1; TRAPPC2
Species	Human
GeneID (Human)	6399
GeneID	6399
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human TRAPPC2 (P0Di81)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P0Di81
Protein Weight	16 kDa
Shipping	Ice bag