

Anti-Trappc2 (Sedl) Rabbit Polyclonal Antibody For Wb, Elisa - P0Di81

Item Number: CM0015992

Introduction

Rabbit polyclonal antibody against human TRAPPC2 (SEDL), 16 kDa. Suitable for Western blot and ELISA. Validated in human samples; immunogen covers full-length 1-140. Ideal for vesicular transport and transcription studies.

[Learn More](#)

Parameter	Value
Product Name	TRAPPC2 Rabbit pAb
Remarks/Alias	SEDL; SEDT; MIP2A; TRS20; ZNF547L; hYP38334; TRAPPC2P1; TRAPPC2
Species	Human
GeneID (Human)	6399
GeneID	6399
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-140 of human TRAPPC2 (NP_055378.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P0Di81
Protein Weight	16kDa
Shipping	Ice bag