

Anti-Trappc6A Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O75865

Item Number: CM0020873



Introduction

TRAPPC6A Rabbit Polyclonal Antibody validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Suitable for melanosome biogenesis and vesicular transport research. One-stop sourcing from CamelBio.

[Learn More](#)

Parameter	Value
Product Name	TRAPPC6A Rabbit pAb
Remarks/Alias	TR533; TRAPPC6A
Species	Human
GeneID	79090
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-173 of human TRAPPC6A (NP_077013.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse, Rat
SWISS	O75865
Protein Weight	18kDa
Shipping	Ice bag