

Anti-Trappc6B Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q86Sz2

Item Number: CM0016118



Introduction

Rabbit polyclonal antibody targeting human TRAPPC6B, suitable for WB, IHC-P, and ELISA. Shows cross-reactivity with mouse and rat. Supports research on intracellular trafficking and early neural circuitry development.

[Learn More](#)

Parameter	Value
Product Name	TRAPPC6B Rabbit pAb
Remarks/Alias	TPC6, NEDMEBA, TRAPPC6B
Species	Human
GeneID (Human)	122553
GeneID	122553
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-130 of human TRAPPC6B (NP_803235.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q86SZ2
Protein Weight	18kDa
Shipping	Ice bag