

Anti-Trem1 Rabbit Polyclonal Antibody For If/Icc, Elisa - Q9Np99

Item Number: CM0019328



Introduction

High-quality rabbit polyclonal antibody targeting human TREM1 (CD354). Validated for IF/ICC and ELISA applications. 26 kDa protein. Ideal for innate immunity and inflammatory response research, including septic shock and atherosclerosis studies.

[Learn More](#)

Parameter	Value
Product Name	TREM1 Rabbit pAb
Remarks/Alias	CD354; TREM-1; TREM1
Species	Human
GeneID (human)	54210
GeneID	54210
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human human TREM1(NP_061113.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q9NP99
Protein Weight	26kDa
Shipping	Ice bag