

Anti-Trex1 Rabbit Monoclonal Antibody For Wb, Ip, Elisa - Q9Nsu2

Item Number: CM0004784



Introduction

Rabbit monoclonal TREX1 antibody validated for WB, IP, ELISA in human. Targets a 33 kDa DNA exonuclease involved in autoimmunity and micronuclei DNA degradation. Essential tool for innate immunity research.

[Learn More](#)

Parameter	Value
Product Name	TREX1 Rabbit mAb
Remarks/Alias	CRV; AGS1; DRN3; HERNS; RVCLS; TREX1
Species	Human
GeneID (Human)	11277
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 150-314 of human TREX1 (Q9NSU2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human
SWISS	Q9NSU2
Protein Weight	33kDa
Shipping	Ice bag