

Anti-Trex1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Nsu2

Item Number: CM0028136



Introduction

TREX1 rabbit polyclonal antibody validated for WB, IF/ICC, and ELISA applications. Detects TREX1 in human, mouse, and rat samples. Targets the 3'-5' exonuclease TREX1 (Q9NSU2), critical for DNA degradation and innate immune regulation.

[Learn More](#)

Parameter	Value
Product Name	TREX1 Rabbit pAb
Remarks/Alias	CRV; AGS1; DRN3; HERNS; RVCLS; TREX1
Species	Human
GeneID (Human)	11277
GeneID	11277
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human TREX1 (NP_057465.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NSU2
Protein Weight	33kDa
Shipping	Ice bag