

Anti-Trh Polyclonal Antibody For Wb, If-P, Elisa - P20396

Item Number: CM0026053

Introduction

Rabbit polyclonal antibody targeting human pro-thyrotropin-releasing hormone (TRH). Validated for WB, IF-P, ELISA. Cross-reacts with human, mouse, rat. Predicted MW ~27 kDa. For research on HPT axis and hair growth.

[Learn More](#)

Parameter	Specification
Product Name	TRH Rabbit pAb
Remarks/Alias	TRF; Pro-TRH; TRH
Species	Human
Gene ID (Human)	7200
Gene ID	7200
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 25-230 of human TRH (NP_009048.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P20396
Protein Weight	27kDa
Shipping	Ice bag