

Anti-Trh Polyclonal Antibody For Wb, If/Icc, Elisa - P20396

Item Number: CM0015326

Introduction

High-quality rabbit polyclonal antibody against human TRH (pro-thyrotropin-releasing hormone). Suited for WB, IF/ICC, ELISA; cross-reacts with mouse and rat. Supports studies in TSH secretion and hair follicle biology.

[Learn More](#)

Parameter	Value
Product Name	TRH Rabbit pAb
Remarks/Alias	TRF; Pro-TRH; TRH
Species	Human
GeneID (human)	7200
GeneID	7200
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 143-242 of human TRH (NP_009048.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P20396
Protein Weight	27kDa
Shipping	Ice bag