

Anti-Trhr Rabbit Polyclonal Antibody For Wb, Elisa - P34981

Item Number: CM0020252



Introduction

High-quality anti-TRHR rabbit polyclonal antibody validated for Western blot and ELISA. Targets human thyrotropin-releasing hormone receptor (P34981), a GPCR involved in thyroid regulation. Ideal for IVD and research applications.

[Learn More](#)

Parameter	Value
Product Name	TRHR Rabbit pAb
Remarks/Alias	CHNG7; TRH-R; TRHR
Species	Human
GeneID (Human)	7201
Immunogen	Recombinant Protein: Recombinant fusion protein containing a sequence corresponding to amino acids 320-398 of human TRHR (NP_003292.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P34981
Protein Weight	45kDa
Shipping	Ice bag