

Anti-Triap1 Polyclonal Antibody For Wb, Elisa - 043715

Item Number: CM0027343



Introduction

Rabbit polyclonal TRIAP1 antibody validated for WB and ELISA, reactive to human and mouse. Targets the 9 kDa mitochondrial apoptosis regulator (p53CSV/WF-1). Ideal for cell death studies.

[Learn More](#)

Parameter	Value
Product Name	TRIAP1 Rabbit pAb
Remarks/Alias	WF-1; MDM35; P53CSV; HSPC132; TRIAP1
Species	Human
GeneID (Human)	51499
GeneID	51499
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-76 of human TRIAP1 (NP_057483.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O43715
Protein Weight	9kDa
Shipping	Ice bag