


Anti-Trib1 Rabbit Polyclonal Antibody For Wb, Elisa - Q96Ru8

Item Number: CM0014166



Introduction

High-specificity rabbit polyclonal antibody targeting human TRIB1 for Western blot and ELISA applications. Useful for myeloid differentiation and CEBPA degradation studies. Detects TRIB1, also known as GIG-2/SKIP1, with a molecular weight of 41 kDa.

[Learn More](#)

Parameter	Value
Product Name	TRIB1 Rabbit pAb
Remarks/Alias	C8FW; GIG2; TRB1; GIG-2; SKIP1; TRB-1; TRIB1
Species	Human
GeneID (Human)	10221
GeneID	10221
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 203-372 of human TRIB1 (NP_079471.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q96RU8
Protein Weight	41kDa
Shipping	Ice bag