


Anti-Trib2 Polyclonal Antibody For Wb, Elisa - Q92519

Item Number: CM0015440



Introduction

Rabbit polyclonal antibody against human TRIB2, validated for western blot and ELISA. Suitable for human and mouse samples; targets the Tribbles homolog 2 pseudokinase involved in MAP kinase signaling regulation.

[Learn More](#)

Parameter	Value
Product Name	TRIB2 Rabbit pAb
Remarks/Alias	C5FW; TRB2; GS3955; TRIB2
Species	Human
GeneID (Human)	28951
GeneID	28951
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of human TRIB2 (NP_067675.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q92519
Protein Weight	39kDa
Shipping	Ice bag