

## Anti-Trib3 Rabbit Monoclonal Antibody For Wb And Elisa - Q96Ru7

Item Number: CM0012081

### Introduction

High-quality rabbit monoclonal antibody targeting human TRIB3. Validated for Western blot and ELISA applications. Ideal for apoptosis and stress response studies.

[Learn More](#)

| Parameter        | Value                                                                        |
|------------------|------------------------------------------------------------------------------|
| Product Name     | TRIB3 Rabbit mAb                                                             |
| Remarks/Alias    | NIPK; SINK; TRB3; SKIP3; C20orf97; TRIB3                                     |
| Species          | Human                                                                        |
| GeneID (Human)   | 57761                                                                        |
| GeneID           | 57761                                                                        |
| Immunogen        | Synthetic peptide corresponding to amino acids 1-100 of human TRIB3 (Q96RU7) |
| Source           | Rabbit                                                                       |
| Category         | Monoclonal Antibodies                                                        |
| Application      | WB, ELISA                                                                    |
| Cross Reactivity | Human                                                                        |
| SWISS            | Q96RU7                                                                       |
| Protein Weight   | 40kDa                                                                        |
| Shipping         | Ice bag                                                                      |