

Anti-Trif / Ticam1 Rabbit Polyclonal Antibody For Wb, Elisa - Q8Iuc6

Item Number: CM0016146



Introduction

Rabbit polyclonal antibody targeting human TRIF/TICAM1 (aa 180-280). Validated for WB and ELISA with cross-reactivity to mouse and rat. Ideal for innate immunity, TLR signaling, and antiviral research.

[Learn More](#)

Parameter	Value
Product Name	TRIF / TICAM1 Rabbit pAb
Remarks/Alias	TRIF; IIAE6; MyD88-3; PRVTIRB; TICAM-1; TRIF / TICAM1
Species	Human
GeneID (human)	148022
GeneID	148022
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 180-280 of human TRIF / TICAM1 (NP_891549.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8IUC6
Protein Weight	76kDa
Shipping	Ice bag