

Anti-Trim17 Polyclonal Antibody For Western Blot And Elisa - Q9Y577

Item Number: CM0020302

Introduction

Rabbit polyclonal antibody targeting human TRIM17 (E3 ubiquitin ligase). Suitable for Western blot and ELISA, cross-reacts with mouse and rat. Ideal for studying ubiquitination, neuronal apoptosis, and autophagy pathways.

[Learn More](#)

Parameter	Value
Product Name	TRIM17 Rabbit pAb
Remarks/Alias	RBCC; terf; RNF16; TRIM17
Species	Human
Gene ID (Human)	51127
Gene ID	51127
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 280-477 of human TRIM17 (NP_057186.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y577
Protein Weight	54kDa
Shipping	Ice bag