


Anti-Trim21/Ss-A Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P19474

Item Number: CM0017799



Introduction

Rabbit polyclonal antibody against human TRIM21/SS-A (RO52, SSA1, RNF81). Validated for WB, IHC-P, IF/ICC, IP, ELISA with cross-reactivity to mouse and rat. 54kDa. Ideal for autoimmunity and ubiquitination research.

[Learn More](#)

Parameter	Value
Product Name	TRIM21/SS-A Rabbit pAb
Remarks/Alias	SSA; RO52; SSA1; RNF81; Ro/SSA; TRIM21/SS-A
Species	Human
GeneID	6737
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P19474
Protein Weight	54kDa
Shipping	Ice bag