


Anti-Trim24 Polyclonal Antibody For Wb, Elisa - O15164

Item Number: CM0016610



Introduction

High-quality anti-TRIM24 rabbit polyclonal antibody validated for WB and ELISA. Targets human TRIM24, a transcriptional coactivator and E3 ubiquitin ligase with key roles in DNA damage response, innate immunity, and cancer. Ideal for IVD developers.

[Learn More](#)

Parameter	Value
Product Name	TRIM24 Rabbit pAb
Remarks/Alias	PTC6; TF1A; TIF1; RNF82; TIF1A; hTIF1; TIF1ALPHA; TRIM24
Species	Human
Gene ID (Human)	8805
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 600-820 of human TRIM24 (NP_056989.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O15164
Protein Weight	117kDa
Shipping	Ice bag