


Anti-Trim24 Rabbit Monoclonal Antibody For Wb, Elisa - O15164

Item Number: CM0003170



Introduction

Highly specific rabbit monoclonal antibody against TRIM24, validated for Western blotting (WB) and ELISA. Detects human, mouse, and rat TRIM24 (117kDa). Ideal for transcription regulation, DNA damage response, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	TRIM24 Rabbit mAb
Remarks/Alias	PTC6; TF1A; TIF1; RNF82; TIF1A; hTIF1; TIF1ALPHA; TRIM24
Species	Human
Gene ID (Human)	8805
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 400-500 of human TRIM24 (O15164).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15164
Protein Weight	117kDa
Shipping	Ice bag