


## Anti-Trim25 Monoclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - Q14258

Item Number: CM0010782



### Introduction

High-specificity rabbit monoclonal antibody targeting human TRIM25 (Q14258), validated for western blot, IHC, IF, IP, and ELISA. Ideal for innate immunity, antiviral signaling, and ubiquitination research. Reliable detection of the 71 kDa protein.

[Learn More](#)

| Parameter        | Value                                                                                                               |
|------------------|---------------------------------------------------------------------------------------------------------------------|
| Product Name     | [KD Validated] TRIM25 Rabbit mAb                                                                                    |
| Remarks/Alias    | EFP; Z147; RNF147; ZNF147                                                                                           |
| Species          | Human                                                                                                               |
| GeneID (human)   | 7706                                                                                                                |
| GeneID           | 7706                                                                                                                |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 190-379 of human TRIM25 (NP_005073.2) |
| Source           | Rabbit                                                                                                              |
| Category         | Monoclonal Antibodies                                                                                               |
| Application      | WB, IHC-P, IF/ICC, IP, ELISA                                                                                        |
| Cross Reactivity | Human                                                                                                               |
| SWISS            | Q14258                                                                                                              |
| Protein Weight   | 71kDa                                                                                                               |
| Shipping         | Ice bag                                                                                                             |