

Anti-Trim25 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - Q14258

Item Number: CM0016903



Introduction

Rabbit polyclonal antibody targeting human TRIM25 (Q14258), an E3 ubiquitin/ISG15 ligase essential for innate immunity. Validated for WB, IHC-P, IF/ICC, IP, ELISA. Cross-reacts with mouse and rat. Ideal for antiviral signaling and cancer research.

[Learn More](#)

Parameter	Value
Product Name	TRIM25 Rabbit pAb
Remarks/Alias	EBP; Z147; RNF147; ZNF147; TRIM25
Species	Human
Gene ID (Human)	7706
Gene ID	7706
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 190-379 of human TRIM25 (NP_005073.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14258
Protein Weight	71kDa
Shipping	Ice bag