

Anti-Trim29 Polyclonal Antibody For Wb, Elisa - Q14134

Item Number: CM0025008



Introduction

Rabbit polyclonal antibody against human TRIM29 (ATDC), validated for Western blot (WB) and ELISA. Reacts with human, mouse, and rat samples. Suitable for innate immunity and macrophage activation research.

[Learn More](#)

Parameter	Value
Product Name	TRIM29 Rabbit pAb
Remarks/Alias	ATDC; TRIM29
Species	Human
GeneID (Human)	23650
GeneID	23650
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 309-588 of human TRIM29 (NP_036233.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14134
Protein Weight	66kDa
Shipping	Ice bag