


Anti-Trim31 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Bzy9

Item Number: CM0014909



Introduction

Rabbit polyclonal antibody targeting TRIM31, an E3 ubiquitin ligase regulating antiviral immune response and inflammation. Validated for WB, IF/ICC, ELISA with human and mouse cross-reactivity. Ideal for innate immunity research.

[Learn More](#)

Parameter	Value
Product Name	TRIM31 Rabbit pAb
Remarks/Alias	RNF; HCG1; HCGI; C6orf13; TRIM31
Species	Human
Gene ID (Human)	11074
Gene ID	11074
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 186-425 of human TRIM31 (NP_008959.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9BZY9
Protein Weight	48kDa
Shipping	Ice bag