

Anti-Trim32 Rabbit Polyclonal Antibody For Wb, Elisa - Q13049

Item Number: CM0026044



Introduction

Rabbit polyclonal antibody against TRIM32 (E3 ubiquitin ligase) for WB, ELISA. Validated in human, mouse, rat. Ideal for studying innate immunity, autophagy, and muscular dystrophy.

[Learn More](#)

Parameter	Specification
Product Name	TRIM32 Rabbit pAb
Remarks/Alias	HT2A; BBS11; TATIP; LGMD2H; LGMDR8; TRIM32
Species	Human
Gene ID (Human)	22954
Gene ID	22954
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human TRIM32 (NP_036342.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13049
Protein Weight	72kDa
Shipping	Ice bag