

Anti-Trim38 Rabbit Polyclonal Antibody For Wb And Elisa - 000635

Item Number: CM0028646

Introduction

Rabbit polyclonal antibody targeting human TRIM38, suitable for Western blot and ELISA applications. Effective for studying innate immune responses and ubiquitin-proteasome pathways. Immunogen corresponds to amino acids 266-465 of human TRIM38.

[Learn More](#)

Parameter	Value
Product Name	TRIM38 Rabbit pAb
Remarks/Alias	RNF15; RORET; TRIM38
Species	Human
Gene ID (Human)	10475
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 266-465 of human TRIM38 (NP_006346.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O00635
Protein Weight	53kDa
Shipping	Ice bag